CHEMISTRY STUDY MATERIALS FOR CLASS 12 (Investigatory Project)-3RD PART GANESH KUMAR DATE:- 23/01/2021

TEST FOR CARBON DIOXIDE

EXPERIMENT: -

As soon as the bottles were opened, one by one the sample was passed through lime

water. The lime water turned milky.

OBSERVATION:

Serial	Name of the	Time Taken	Conclusion
Number	Drink	(in second)	
01	COCA COLA	26.5	CO_2 is present
02	SPRITE	21	CO_2 is present
03	LIMCA	35	CO_2 is present
04	FANTA	36	CO_2 is present
INFEREN	ICE:		

All the soft drinks contain dissolved Carbon dioxide in water. The Carbon dioxide dissolves in water to form Carbonic acid, which is responsible for its tangy taste.



EXPERIMENT: -

Glucose is a reducing Sugar acid. Its presence is detected by the following test:

1. BENEDICT'S REAGENT TEST: -

Small samples of cold drinks of different brands were taken in a test tube and a few drops of Benedict's reagent were added. The test tube was heated for few seconds . Formation of reddish colour confirmed the presence of glucose in cold drinks.

OBSERVATION:

Serial Number	Name of the Drink	Observation	Conclusion
01	COCA COLA	REDDISH COLOUR	GLUCOSE IS PRESENT
-		PRICIPETATE	
07	SPRITE	REDDISH COLOUR	GLUCOSE IS
U Z		PRICIPETATE	PRESENT
03	LIMCA	REDDISH COLOUR	GLUCOSE IS
05		PRICIPETATE	PRESENT
01	FANTA	REDDISH COLOUR	GLUCOSE IS
04		PRICIPETATE	PRESENT

INFERENCE: All the samples gave positive test for Glucose with

Benedict's reagent. Hence all the drinks contain Glucose.

2. FEHLING'S SOLUTION TEST

EXPERIMENT: -

Small samples of cold drinks of different brands were taken in a test tube and a few drops of Fehling's reagent A solution and Fehling's B solution was added in equal amount. The test tube was heated in a water bath for 10minutes. Appearance of brown precipitate confirmed the presence of Glucose in Cold drinks.

OBSERVATION:

Serial Number	Name of the Drink	Observation	Conclusion
01	COCA COLA	REDDISH BROWN	GLUCOSE IS
		PRECIPITATE	PRESENT
02	SPRITE	REDDISH BROWN	GLUCOSE IS
		PRECIPITATE	PRESENT
03	LIMCA	REDDISH BROWN	GLUCOSE IS
		PRECIPITATE	PRESENT
04	FANTA	REDDISH BROWN	GLUCOSE IS
• •		PRECIPITATE	PRESENT

INFERENCE:

All the samples gave a positive test for Glucose with Fehling's (A & B) solutions. Hence all the cold drink contains Glucose.

<u>TEST FOR PHOSPHATE</u>

EXPERIMENT: -

Small samples of each brand of cold drinks were taken in separate test tubes and Ammonium Molybdate followed by concentrated Nitric Acid was added to it. The solutions were heated. Appearance of canary-yellow precipitate confirmed the presence of Phosphate ions in Cold drinks.

OBSERVATION:			
Serial Number	Name of the Drink	Observation	Conclusion
01	COCA	CANARY-YELLOW	PHOSPHATE IS
	COLA	PRECIPITATE	PRESENT
02	SPRITE	CANARY-YELLOW	PHOSPHATE IS
		PRECIPITATE	PRESENT
03	LIMCA	CANARY-YELLOW	PHOSPHATE IS
		PRECIPITATE	PRESENT
04	FANTA	CANARY-YELLOW	PHOSPHATE IS
		PRECIPITATE	PRESENT

INFERENCE: All the samples gave positive test for Phosphate.

Hence all the Cold drinks contain Phosphate ions.

CHEMICALREACTION INVOLVED TEST FOR & LCOHOL

EXPERIMENT: -

Small samples of each brand of cold drinks were taken in separate test tubes and Iodine followed by Potassium Iodide and Sodium Hydroxide solution was added to each test tube. Then the test tubes solutions were heated in water bath for 30 minutes. Appearance ofyellow precipitate confirmed the presence of Alcohol in Cold drinks.

Serial Number	Name of the Drink	Observation	Conclusion
01	COCA COLA	YELLOW PRECIPITATE	ALCOHOL IS PRESENT
02	SPRITE	YELLOW PRECIPITATE	ALCOHOL IS PRESENT
03	LIMCA	YELLOW PRECIPITATE	ALCOHOL IS PRESENT
04	FANTA	YELLOW PRECIPITATE	ALCOHOL IS PRESENT

OBSERVATION:

INFERENCE: All the samples gave positive test for Alcohol. Hence all

the Cold drinks contain Glucose. CH₃CH₂OH+4I₂+6NaOH+CHI₃+HCOONa+5NaI+
